



RIDS panel

For the screening of susceptibility to respiratory infectious diseases

RIDS panel

RIDS panel (Respiratory Infectious Diseases Susceptibility) is a kit designed for large population screening through a molecular protocol based on NGS technologies. The resulting genetic profile evaluates both genetic variants in the infection pathway (predisposition or protection against infection), and variants of the genes involved in the immune response to infection (predisposition or protection against a worse outcome).

The kit is validated for analysis of DNA extracted from different body tissues (blood, saliva, etc.).

TECHNOLOGY

RIDS panel kit contains all reagents required for the preparation of a specific bidirectional library of amplicons designed for the NGS analysis using Ion Torrent sequencers.

WORKFLOW

The RIDS panel kit is part of a DNA-to-variant solution that offers streamlined content, easy-to-perform library preparation, push-button sequencing systems, and simplified data analysis.

Library preparation follows a straightforward, PCR-based protocol that can be completed in as little as 5 hours, with < 1 hours hands-on time. Resulting libraries can be normalized, pooled, and then loaded for sequencing.

SOFTWARE ANALYSIS

4eVAR is our proprietary analysis software. The analysis is designed based on the characteristics and technologies of the kit, in order to increase accuracy of results, and to have the complete control on the entire process.

SAMPLE PER RUN

Instrument	Samples per run*
	Germline
Ion 316™ Chip/Ion 510™ Chip	36
Ion 318™ Chip/Ion 520™ Chip	80
Ion 530™ Chip	240
Ion PI™ Chip/Ion 540™ Chip	720

*the estimated maximum number of samples per chip assumes a reading depth of 300x for the germline. The optimal number of samples can be empirically estimated on the local setup.

ORDERING INFORMATION

Product**	REF
RIDS panel	H1080-16 (16 test)
RIDS panel	H1080-96 (96 test)
Barcode	
Barcode Set 1-16	6001
Barcode Set 17-32	6002
Barcode Set 33-48	R6003
Barcode Set 49-64	R6004

**the kit is also available in its version only for research use (RUO).

Table 1: List of hotspots in RIDS panel

		gene	RSID
			ABO (depends on group)
susceptibility/resistance	aspecific	DDP4 (associated with resistance)	rs13015258 rs117888248 rs116302758 rs56179129 rs115450134
		CD147 (BSG)	rs201850688 rs11551906 rs144824657 rs41276870
		CCR5 (associated with resistance)	rs333
		ACE1	rs4341
	GW specific for COVID19	ACE2 (associated with resistance to infection)	K31R, N33I, H34R, E35K, E37K, D38V, Y50F, N51S, M62V, K68E, F72V, Y83H, G326E, G352V, D355N, Q388L, and D509Y
		TMPRSS2 (associated with increased susceptibility)	rs2070788 rs383510 rs200291871 rs75603675 rs61735791 rs114363287 rs12329760
		GLL5, GNAZ, RSPH14, RAB36 and BCR	rs73166864
		IVNS1ABP, SWT1	rs6668622
		ApoE	rs429358-C-C (e4e4)
		IFITM3	rs12252 rs6598045
outcome (associated with worst prognosis)	SLC6A20, LZTFL1, CCR9, FYCO1, CXCR6, XCR1	rs11385942 rs73064425	
	OAS3	rs10735079	
	TMEM189-UBE2V1	rs6020298-A	
	DPP9	rs2109069	
	PCSK3 (associate with)	rs16944971 rs780909157 rs201551785 rs769208985 rs1236237792	
	SRRM1,IVNS1ABP (hospitalization)	rs111972040	